

MEDI LIFTEEZ

NOW ON FORMULARY

MEDI LIFTEEZ Adhesive Remover Non-Sting

Benefits:

- Effectively remove dressings, pouches and medical devices, minimising the potential for skin stripping and pain on dressing removal
- Fast action - starts breaking adhesive bonds on application
- Silicone based - Alcohol free
- Non-Sting formulation, pain free application
- Dries in seconds
- Aerosol can be sprayed and used at any angle

**AEROSOL
& WIPE
FORMATS**



Types of MARS I		Appliances LIFTEEZ can be used with
	<p>Vesicle Localised epidermal elevations (tiny blisters) in the skin containing clear fluid and less than ½ cm in diameter</p>	<ul style="list-style-type: none"> • Small or large adhesive dressings • NPWT - Negative Pressure Wound Therapy • Medical/Surgical tapes • Adhesives securing devices such as central lines, NG tubes, IV and SC cannula • Ostomy pouches and stoma appliance • Urinary sheaths • ECG electrodes • Theatre adhesive drapes <p><i>Note: MARS I - Medical Adhesive Related Skin Injury or Skin Stripping</i></p>
	<p>Bulla Localised epidermal elevations (tiny blisters) in the skin containing clear fluid and greater than ½ cm</p>	
	<p>Erosion Erosion is very superficial loss of some or all of the epidermis (the outer layer) leaving a denuded surface (i.e. a derroofed blister)</p>	
	<p>Tear Separation of the epidermis from the dermis, or separation of both the epidermis and dermis from underlying structures</p>	
	<p>Skin Stripping Unintended damage to and loss of the upper layer of the epidermis</p>	

Ordering Information:

Format	Selling Unit	Product Code	PIP Code	NHS-SC Product Code
Wipe	30 Pack	60321	341-3309	EXC042
Aerosol - 50ml	1	60833	389-7147	EXC041

If you require more information about **MEDI LIFTEEZ** or would like a representative to contact you, please email us at info@medicareplus.co.uk or visit us at www.medicareplus.co.uk

